



EL-SMS-2G-VAC

Vaccine Temperature Monitor with SMS Alerts

Receive SMS alerts as soon as the temperature triggers an alarm.

- High and low temperature alarm level notifications
- Daily, weekly or monthly summary SMS messages
- Setup using the free EasyLog SMS Apps available from Google Play and the App Store
- -40 to 60°C (-40 to 140°F) measurement range using supplied probe
- Glycol bottle probe mimics vaccine in storage for accurate temperature monitoring.
- Powered by four Alkaline AA batteries



EasyLog SMS Mobile Apps available from Google Play and the App Store


The EasyLog SMS App can be downloaded for free from Google Play or the App Store, allowing you to manage your EL-SMS devices with ease.

- Animated walkthrough of the setup process
- Easily configure your EL-SMS unit and set temperature alarms
- Receive alarm notifications & regular summary messages
- Historical view of all messages received from your devices

Within the app you can set the temperature alarm levels, set the frequency of summary messages (which contain minimum, maximum, and average readings since the last summary), as well as defining up to three contacts to receive SMS notifications.

Once the options are set the App will guide you through the process of sending a configuration message to the device. Your EL-SMS-2G-VAC is now monitoring - ready to send alarm and summary notifications.

Alarm and summary messages can be sent to any device capable of receiving SMS messages, but the EasyLog SMS App enables you to monitor multiple EL-SMS-2G-VAC devices in one easy package.



SPECIFICATIONS

Glycol Probe <i>(Supplied)</i>	Measurement range	-40 to 60°C (-40 to 140°F)
	Accuracy (-10 to 60°C / 14 to 140°F)	±0.1°C (±0.2°F)
	Accuracy (-40 to 60°C / -40 to 140°F)	±0.3°C (±0.6°F)
Accuracy (Monitor error)		±0.1°C (±0.2°F)
Internal resolution		0.1°C (0.1°F)
Operating temperature range		-20 to 40°C (-4 to 104°F)
Battery life		1 year (at 25°C and daily summary message)
Dimensions		105 x 105 x 29mm (4.1 x 4.1 x 1.1")
Frequency Band		Quad-Band 2G GSM (850/900/1800/1900MHz)

ACCESSORIES

BAT 1V5 AA	Replacement battery (4 required)
EL-P-VAC	Replacement glycol thermistor probe
EL-PROBE-EXTENDER-5.0M-TP	Probe extension (5.0m)
EL-PROBE-EXTENDER-10.0M-TP	Probe extension (10.0m)

INCLUDED IN THE BOX

x4 BAT 1V5 AA	Batteries
EL-P-VAC	Glycol thermistor probe
M-EL-SMS-WM	Mounting Bracket
CAL-TEMP-UK	Temperature Calibration certificate



CALIBRATION CERTIFICATES INCLUDED

The EL-SMS-2G-VAC comes with a calibration certificate issued at 2 and 8°C, using reference equipment which has been calibrated by a UKAS/NIST accredited laboratory and using apparatus traceable to national or international standards.




EL-SMS-2G-VAC



Vaccine Temperature Monitor with SMS Alerts

Receive SMS alerts as soon as the temperature triggers an alarm.

EL-SMS-2G-PROBE-G



- ① Status button
- ② Stop button
- ③ Status LEDs
- ④ Netlight
- ⑤ External temperature probe port
- ⑥ Sounder
- ⑦ Battery Bay
- ⑧ Reset switch
- ⑨ SIM card holder

SIM CONTRACT

To use your EL-SMS-2G-VAC you will need to install a pre-registered micro SIM card from a 2G enabled service provider of your choice. It is important that no PIN code is enabled as this will prevent the unit from operating.


Visit maps.mobileworldlive.com to check for GSM coverage in your geographic location.

BATTERY INFORMATION

For reliable operation of the unit it is recommended that only high quality alkaline batteries, such as those supplied, are used.

Please refer to the EasyLog SMS User Manual for more information.

THERMISTOR PROBE ACCURACY



Supplied 'Glycol' Thermistor Probe
Part Number: EL-P-VAC



Issue 3 07/2019 Page 2 of 2

CMI BELGIUM

Avenue de Ramelot 3
B - 1480 Saintes (Sint Renelde)
Tel : +32 (0)2 531 11 21
info@cmitest.com

CMI THE NETHERLANDS

Luchthavenweg 81
NL - 5657 EA Eindhoven
Tel : +31 (0)402 88 87 57
info@cmitest.com

CMI FRANCE

Ecoparc d'Affaires
F- 41210 Neung-sur-Beuvron
Tel : +33 (0)2 54 95 70 95
info@cmitest.com

Agence Rhône Alpes :
F- 26160 La Bégude-de-Mazenc
Tel : +33 (0)4 75 54 57 26



www.cmitest.com