

## **Cleanroom Management International** Groupe Trescal







#### CRĒDO™ OPERATING ROOM CONTAINERS

Blood banks, hospitals, emergency medical services, and military organisations know the value of protecting blood supplies. Safely transporting materials can mean the difference between life and death and professionals need a solution that will ensure precise thermal control with the ability to be seamlessly reinstated into inventory.

At the core of every Crēdo™ Operating Room Container is our patented TIC™ system surrounding the critical payload with walls of PCM (phase change material) to maintain required temperature range and duration.

#### WHY CHOOSE CRĒDO™ OPERATING ROOM CONTAINER?

- Patented TIC™ system and PCMbased
- Reusable phase change technology reduces costs and eliminates the complexity and restrictions of outdated components such as wet/dry ice and gel packs
- Designed and tested to last in harsh conditions
- Quality, durable container construction



## **GOLDEN HOUR™ OPERATING ROOM CONTAINER**



**DURATION** 10 Hours



**TEMPERATURE** 

**SERIES 4** 1°C to 6°C **SERIES 22** 20°C to 24°C



্বা VOLUME

## CRĒDO™ OPERATING **ROOM CONTAINER**



**DURATION** 16 Hours



**SERIES 4 - 824 TEMPERATURE** 1°C to 6°C




□ ✓ VOLUME

# CRĒDO™ OPERATING ROOM CONTAINERS

	Payload Size	Litres	<b>Exterior Dimensions</b> WxDxH		Interior Dimensions WxDxH		Tare Weight	
			mm	inchs	mm	inchs	kg	lbs
Series	Series 4 ORC	2L	330 x 241 x 216	13 x 10 x 9	152 x 127 x 108	6×5×4	4.1	9
	Series 22 ORC	2L	330 x 241 x 216	13 x 10 x 9	152 x 127 x 108	6×5×4	4.1	9
Series	Series 4 824 ORC	8L	406 x 260 x 381	16 x 10 x 15	282 x 165 x 165	11 x 7 x 7	6.4	14

Whether transporting platelets, RBCs, FFP, cord blood, as well as other blood products for surgical and EMT purposes, Crēdo™ thermal protection containers offer the safest and most reliable temperature-control transport solution.

### **ADDITIONAL DATA**





2-Litres



8-Litres

### **CMI FRANCE**

Ecoparc d'Affaires F- 41210 Neung-sur-Beuvron Tel : +33 (0)2 54 95 70 95 Agence Rhône Alpes : F- 26160 La Bégude de Mazenc Tel: +33 (0)4 75 54 57 26

infofr@cmitest.com www.cmitest.com





